

Certificate of Analysis

GHK Copper 100mg


copper;(2S)-6-amino-2-[[[(2S)-2-[(2-aminoacetyl)amino]-3-(1H-imidazol-5-yl)propanoyl]amino]hexanoate

Compound : GHK-Cu
Lot number : 2025-04-07
Analysis date : 2025-05-08
Purity % : 99.97%
Quantity : 139.82 mg
Method : HPLC-UV-MS

Client : Transforma Peptides
www.transformapeptides.com

PubChem CID: 71587328
<https://pubchem.ncbi.nlm.nih.gov/compound/71587328>

High Performance Liquid Chromatography (HPLC) UV – Quantitative test



PEAK LIST			
Number of detected peaks: 3			
	Time (min)	Area	%Area
1	4.77	5.69E+01	0.03
2	6.23	1.65E+05	99.65
3	6.35	5.25E+02	0.32

GHK Cu
GHK

Overall purity : 99.97

The GHK-Copper complex dissociates slightly during HPLC resulting in the small GHK peak immediately after GHK-Cu. This is normal.

Quantification by HPLC-UV
Measured quantity : 139.82 mg/vial

Note: Injectable peptides may contain salts and sugars to aid in solubility and act as pH buffers. These are not normally detected using UV and are not considered impurities.

Analysis Performed by
Ken Pendarvis, ChE
Analytical Chemist
MZ Biolabs
contact@mzbiolabs.com



2025-05-09

GHK Copper 100mg

PubChem CID: 71587328

<https://pubchem.ncbi.nlm.nih.gov/compound/71587328>

Mass Spectrometry (MS) – Identity Test

Identity confirmed using HPLC-MS

Molecular weight calculated using monoisotopic m/z values from mass spectrum

Expected monoisotopic mass : 402.10 Da


Measured monoisotopic mass : 402.11 Da

Molecular weight confirmed

Note : Monoisotopic m/z values are not easily seen in full spectrum view for larger molecules and peptides.

The dominant isotopic peak (base peak) shown in the spectrum below can be used to approximate the average molecular weight frequently reported by vendors and databases as a secondary means of confirmation.

Recorded MS spectrum



Analysis Performed by
Ken Pendarvis, ChE
Analytical Chemist
MZ Biolabs
contact@mzbiolabs.com



2025-05-09